

# A15 Automated External Defibrillator



# A15

## HeartOn AED

**ANYTIME, ANYWHERE,  
ANYONE**

Mediana's HeartOn AED helps to save people's lives when Sudden Cardiac Arrest occurs.



## CONFIGURATION

### 1 Status Indicator

Unit status, battery level status and temperature range status

### 2 Connecting Socket

Defibrillator pads connector socket  
The pads are always ready for use.

### 3 Visual Icons

Indicator LED lights flash with the corresponding voice prompts. Visual Icons clarify the necessary actions.

### 4 Shock Button

Flashing button indicates "Ready for shock delivery". Push the button to deliver shock.

### 5 Adult/Pediatric Mode Switch

Easily switch from adult to pediatric mode without changing pads.

### 6 Open Switch

Slide the switch to the right to open the cover. Unit automatically turns on.

### 7 Speaker

Acoustic indicator for real time guidance.



IR Port

SD Card Port



Battery

Disposable LiMnO2 Non-Rechargeable



Adult/Pediatric Pads

2-in-1 Electrode


# FEATURES & FUNCTIONS




## 1. Reduction of the defibrillation time

- In order to shorten the time before administering the shock, ECG analysis & charging are performed simultaneously.

### »» Before



### »» After



## 2. Auto-Start Function

- The device is activated by opening the device cover. No need to press a power button.

## 3. Fast response with 2-in-1 Pads

- Adult or Pediatric modes can be selected by simply changing a switch.
- Emergency response time can be shortened, as the right pad is connected to the device at all times.




## 4. Intuitive

- Anyone can use the device with the guidance of icons and voice prompt.


## HOW TO USE




Upon opening the cover, the power is on (Auto-Start)



Attach the pads






Simultaneous ECG analysis & charging



Push the shock button

## AED EVENT REVIEW SOFTWARE



### Mediana Seoul Office

53 Donggwang-ro, Seocho-gu, Seoul,  
South Korea, 06561

Tel : +82 2 542 3375 Fax : +82 2 542 7447

### Mediana Head Office / Factory / R&D

132 Donghwagongdan-ro, Munmak-eup, Wonju-si,  
Gangwon-do, South Korea, 26365

Tel : +82 33 742 5400 Fax : +82 33 742 5483

[www.mediana.co.kr/eng](http://www.mediana.co.kr/eng)

